



# Programming your Minuet Watch

Green cancel key



**Step 1.** Press and hold the **green cancel key** until the unit beeps.

The red light will start flashing.

Release the green cancel key. **Step 2.**

Repeat Step 1

**Step 3.** Press and release the **alarm button** on your new pendant.

The unit will beep to confirm acceptance.

**Step 4.** Press and release the **green cancel key**.

The unit will beep.

The red light should stop flashing.

(If the red light continues to flash, press the green cancel key)

**Step 5.** Make a **test call using** the new programmed **Minuet Watch**.



Actual size

If you have any problems please contact us on 1850 23 23 24

**EMERGENCY  
RESMNSE**

All the reassurance you need

**Tunstall**